American Made Quality

World Class Customer Service



Call Now For Your Free Professional Evaluation! 800-475-8122



A Word About Our Policies...

WARRANTY: This Limited Warranty is extended to the original purchaser/user and cannot be transferred. All Transfer Master Hi-Low beds and accessories are warranted for a period of One Year against defects. A detailed version of this warranty accompanies each bed and is available upon request.

Warranty does not cover products damaged in shipping. It is the responsibility of the customer to inspect the delivered goods upon arrival. Should damages exist, report it immediately to Transfer Master Products, Inc. and the carrier. The carrier will be responsible for repairing or replacing the product. Failure to report these damages upon delivery could result in a loss to the customer.

SHIPPING: All shipments are FOB origin unless otherwise approved by management. All shipments will be shipped using the most economical method unless purchaser requests otherwise and absorbs additional freight charges.

FREIGHT DAMAGE: We package all items very carefully to minimize the likelihood of freight damage. The carrier who delivers our merchandise to you is responsible for loss and damage. Acceptance of the shipment by you is acknowledgement that all of the items were delivered in good condition. Therefore, all claims for loss or damage must be filed immediately with the "FREIGHT CARRIER." Notify our company, as well, so that we may document our records accordingly. When receiving a shipment, always note any visible damage on the delivery ticket. Should you need assistance with any claim, call Customer Service at: 563-864-7364.

CREDIT/RETURNS: All returns must be pre-approved in advance by Transfer Master Products, Inc. and a Return Authorization Form must be returned with product to insure proper credit. Customer is responsible for freight when returning product. Restock Fee is 35 percent plus all freight costs for standard items. Custom, used or altered items may not be returned. Any mattress or upholstered products may not be returned. Transfer Master Products, Inc. reserves the right to refuse any returns after 30 days.

SPECIAL ISSUES

- Dimensions and heights/adjustments stated are approximate. Actual dimensions may vary slightly
- Fabrics, quilt patterns, hand controls, bed styles, etc. may not be shown in photos due to changes beyond our control
- Transfer Master beds are ECU compatible (Environmental Control Unit compatible)
- Due to the custom nature of the products, orders may not be canceled once processing has begun
- If you do not have carpeted floor, ask about our optional glides or caster cups
- Headboard and footboard cannot be used with a dual king if beds are intended to be placed side by side
- Headboard and footboards are not to be used as weight bearing or for repositioning
- Most Hoyer type lifts are compatible with our beds. We cannot guarantee that your Hoyer will work with our beds. TO BE SURE, inquire at the time of placing an order. We may ask for the clearance and width needed for your specific lift
- We manufacture our beds to your requested specifications and within factory tolerances. We cannot guarantee that your specifications
 will work for you. If you have any doubts about your requirements, please consult with a physical therapist or specialist before
 ordering our products
- For your convenience, all of our beds include a 12 foot hand control cord, 2 side retainer bars, and a foot retainer bar
- We reserve the right to refuse service when the product is found to be in an unsanitary condition, soiled, stained or when failure is due to causes other than defective materials or workmanship.

MATTRESS WARRANTY

We offer a wide range of mattresses from cloth orthopedic to computerized turn systems. There are different Warranties depending on the product and often these Warranties are supported by our different providers. Specific Warranties are available upon request. All guarantees are limited to workmanship defects and exclude comfort and other subjective aspects. Customer will be responsible for any freight if the mattress needs to be transported. See policies above for other conditions regarding mattresses. If you purchased your bed from a local Dealer please contact to evaluate any defect.

Innerspring Coiled Cloth Mattresses

Transfer Master beds come standard with an innerspring coiled cloth mattress. As we can not guarantee comfort for any individual, we have weight categories that we use to guide which level of firmness may best suite you:

- Easy Sleep-Approximately 8" thick. Recommended for users under 150 pound
- Supreme Sleep-Approximately 9" thick. Recommended for users 155 pounds to 220 pounds
- Max Sleep-Approximately 10" thick. Recommended for users 225 pounds and higher

Mattress Specifications

- Bonnell coil system
- Motion reducing foam
- Quilted polyester/cotton blend
- Crush resistant edges
- Foam Encased



Quick Facts

- ❖Super Low 9.5 inches from deck to the floor
- ❖Height range 12 inches
- ❖300 pound weight capacity
- Robust, modular control box design
- ❖Tool-free, easy setup split bed
- Mattress included
- ♦Single or dual hand controls
- ◆ECU compatible**
- ❖Compatible with most Hoyer Lift systems**

Available Sizes**

Twin 38"x80" Twin SL 38"x84" 53"x80" Full Full SL 53"x84" 60"x80" Queen Also available in Dual King

Standard 2 or 3 Function Package

- ❖Head Adjustable
- ❖Foot Adjustable
- ❖Hi-Low Adjustable (2 Function model will not have this operation)
- ♦6 Button Hand Control



Standard 4 Function Package

- ❖Head Adjustable
- ❖Foot Adjustable
- ❖Hi-Low Adjustable
- ❖Reverse Trendelenburg
- ❖8 Button Sleek Hand Control



Deluxe 3 Function Package

- ❖Head Adjustable
- ❖Foot Adjustable
- ❖Hi-Low Adjustable
- ♦6 Button Sleek Hand Control
- ❖Deluxe Illuminated Hand Control
- ❖Single Variable Massage (operated only by Deluxe Control)
- ❖Soft Illuminated Floor Light Feature





NIGHT RIDER



Standard 2 or 3 Function Package

- ❖Head Adjustable
- ❖Foot Adjustable
- ❖Hi-Low Adjustable
- ♦6 Button Hand Control



Quick Facts

- Less than 12 inches from deck to the floor
- ♦ Height range 12 inches
- ❖300 pound weight capacity
- Robust, modular control box design
- ❖Tool-free, easy setup split bed
- Mattress included
- ❖Single or Dual hand controls
- ◆ECU Compatible**
- ❖Compatible with most Hoyer Lift systems**

Available Sizes**

Twin 38"x80"
Twin SL 38"x84"
Full 53"x80"
Full SL 53"x84"
Queen 60"x80"
Also available in Dual King

Deluxe

3 Function Package

- ❖Head Adjustable
- ❖Foot Adjustable
- ❖Hi-Low Adjustable
- ❖6 Button Sleek Hand Control
- ❖Deluxe Illuminated Hand Control
- Single Variable Massage (operated only by Deluxe Control)
- ❖Soft Illuminated Floor Light Feature









Demonstrating easy split technology for easy delivery and setup

NEW VALIANT



Quick Facts

- Less than 12 inches from deck to the floor
- ♦ Height range- 12 inches
- ❖5 Function Model
- ❖300 pound weight capacity
- Robust, modular control box design
- ❖Tool-free, easy setup split bed
- Mattress included
- ❖Single or Dual hand controls
- ◆ECU Compatible**
- ❖Compatible with most Hoyer Lift systems**

Available Sizes**

Twin 38"x80" Twin SL 38"x84" Full 53"x80" Full SI 53"x84" Queen 60"x80" Also available in Dual King



5 Function Package ❖Head Adjustable

Standard

- ❖Foot Adjustable
- ❖Hi-Low Adjustable
- Trendelenburg
- ❖Reverse Trendelenburg
- ♦8 Button Sleek Hand Control



Deluxe

5 Function Package

- ❖Head Adjustable
- ❖Foot Adjustable
- ❖Hi-Low Adjustable
- Trendelenburg
- ❖Reverse Trendelenburg
- ❖8 Button Sleek Hand Control
- ❖Deluxe Illuminated Hand Control
- ♦Single Variable Massage (operated only by Deluxe Control)
- ❖Soft Illuminated Floor Light Feature







Demonstrating easy split technology for easy delivery and setup

VALIANT HD/SHD



Quick Facts

- Less than 13.5 inches from deck to the floor
- ❖Height range- 10 inches
- **❖**5 Function
- ♦600/750 pound weight capacity
- ❖Robust, modular control box design
- ❖Tool-free, easy setup split bed
- Mattress included
- ♦Single or Dual hand controls
- ◆ECU Compatible**
- ❖Compatible with most Hoyer Lift systems**

Available Sizes**

Twin 38"x80"
Twin SL 38"x84"
Full 53"x80"
Full SL 53"x84"
Queen 60"x80"
Also available in Dual King



- ❖Looks the same as the Valiant HD/SHD
- ❖3 Function Only
- ❖Also Available in Deluxe Package

Standard 5 Function Package

- ♦ Head Adjustable
- ❖Foot Adjustable
- ❖Hi-Low Adjustable
- Trendelenburg
- ❖Reverse Trendelenburg
- ❖8 Button Sleek Hand Control



Deluxe 5 Function Package

- ♦ Head Adjustable
- ❖Foot Adjustable
- ❖Hi-Low Adjustable
- ❖Trendelenburg
- ❖Reverse Trendelenburg
- ❖8 Button Sleek Hand Control
- ❖Deluxe Illuminated Hand Control
- Single Variable Massage (operated only by Deluxe Control)
- ❖Soft Illuminated Floor Light Feature













Demonstrating easy split technology for easy delivery and assembly for HD/SHD beds

HOSPITAL GRADE HD/SHD



Quick Facts

- Less than 13.5 inches from deck to the floor
- ❖2, 3, or 5 Function packages available
- ❖Height range- 10 inches
- ❖600/750 pound weight capacity
- Powder coat finish for easy cleaning
- Hospital Grade power cord
- ❖5" Swivel lock casters
- ❖Robust, modular control box design
- ❖Tool-free, easy setup split bed
- ♦Single or Dual hand controls
- ◆ECU Compatible**
- ❖Compatible with most Hoyer Lift systems**

Available	Sizes**
Twin	38"x80"
Twin SL	38"x84"
Full	53"x80"
Full SL	53"x84"
Queen	60"x80"







OPTIONS AVAILABLE



Deluxe Side Rails Embedded Controls



Optional Egress Floor Light feature



Hospital Grade Power Strip (at foot of bed)



Straight Wall Bumper



See how easy it is to use a Hoyer-type system with our beds! Hoyer model shown is a 4.5" clearance model.

BED TERMINOLOGY



Hi-Low (Fully Electric/3 Function)

Image illustrates our adjustable bed, shown in fully raised position, demonstrating the ability to vertically raise the bed from its lowest position for ease of transfers, and this also eliminates back strain for care givers.

Head & Foot Adjustable (Semi Electric/2 Function)

Image shows head and foot adjusted (and hi-low fully raised). All Transfer Master beds are head and foot adjustable. This allows for your head and feet to be elevated.



Trendelenburg

The patient is in a supine position, lying face up, and the bed is angled so that the feet are above the head. This position has many benefits, typically with issues such as low blood pressure.



Reverse Trendelenburg

The patient is laid flat in the head-up, feet-down position. This facilitates respiration, and alleviates regurgitation problems and acid reflux.



This helps to elevate the patient from the waist up, which provides relief to the lungs and circulation, and assists him or her in coming from a fully flat position to a sitting position without causing harm or undue strain. The cardiac chair position is usually easier to adjust, with separate sections that support back, bottom, arms, and legs.



Dual King

Our Dual King option is the perfect solution to meet the different needs of any couple. Each bed, separately controlled, can look the same but have entirely different functions. When both beds are in the lowest position and covered with linens, the dual king looks like an Eastern King size bed.



Hover Lift/Patient Lift

This is an assistive device that allows patients in hospitals and nursing homes and those receiving home health care to be transferred between a bed and a chair or other similar resting places, using hydraulic power. Sling lifts are used for patients whose mobility is limited.



ECU Compatible

Environmental Control Unit enables third party control systems (ie. Touchscreens or voice operated systems) to be plugged in to our control boards. Two ports for maximum flexibility!

AVAILABLE OPTIONS

FOR ALL MODELS



Battery Backup allows for bed to be repositioned up to 3 short cycles during a power outage





Joystick Hand Control (New Valiant, Valiant HD/SHD, Hospital Grade Valiant) has large bright colored buttons and a joystick to assist patients with limited mobility to operate their bed with ease



Chrome Half Rails Standard brackets for half rails or 7" extension brackets



Big Button Hand Control (New Valiant, Valiant HD/SHD, Hospital Grade Valiant) has large bright colored buttons to assist patients with limited hand and finger mobility to operate their bed with ease



Headboard & Footboard**
Style as shown



Vertical Deluxe Illuminated Hand Control



Horizontal Deluxe Illuminated Hand Control



Catheter Hook



Caster Cups optional for HD/SHD models



Protection Covers are machine washable, stain protectant, waterproof, easy to fit with dust mite and pet dander protection, cotton toweling, breathable. Available in most sizes.



2" Glides optional for HD/SHD models

PressureGuard® Low Air Loss Systems


Easy Air LR

Unmatched Moisture Control – Now with Lateral Rotation

The Easy Air series offers superior low air loss (see sidebar) and alternating pressure. And now, Easy Air LR delivers a powerful, additional modality: Continuous Lateral Rotation Therapy (CLRT). The unique PressureGuard approach to CLRT has been documented to help wounds close more than three times faster than other commonly provided air therapy surfaces¹.


Supports 500 lbs. (350 lbs. in lateral mode).

1- Anderson, C. et al. Lateral rotation mattresses for wound healing.




Science.

Easy Air's patented design (US patent #6,782,574) maximizes removal of excess moisture. Noncollapsible "Air Diffusion Matrix™" components ensure uninterrupted path for air flow and moisture removal.



PressureGuard Easy Air

surpasses all leading competitors in this key



Data Source: A Proposed Method for Quantifying Low-Air-Loss Mattress Performance by Moisture Transport. Richard S. Figliola PhD, PE, Clemson University, Clemson, SCI. Ostomy/Wound Management. 2003;49(1):32-41.

Ostomy/Wound Management April 2004:50(4):50-62.

Item#	Description
L8039LR-29	Easy Air LR System 80"L x 39"W x 7"H
L8053LR-29	Easy Air LR System 80"L x 53"W x 7"H
8200	Easy Air LR Control Unit

Bariatric Models

Easy Air® XL

Bariatric Low Air Loss for Users up to 1000 lbs.

Heavy Duty low air loss with alternating pressure. Minimizes excess skin moisture – especially important for at-risk obese patients.

ltem#	Description
L8042XL-29	Easy Air XL system, 80"L x 42"W x 7"H
L8048XL-29	Easy Air XL system, 80"L x 48"W x 7"H
* L8053XL-29	Easy Air XL System, 80"L x 53"W x 7"H
L8453XL-29	Easy Air XL system, 84"L x 53"W x 7"H
8300	Easy Air XL Bariatric Control Unit Only

^{*} Lead times: Allow minimum 10 working days from the time of order for shipment.

Geo-Mattress Atlas® with Heel Slope[™] Protection

Bariatric Therapeutic Foam Mattress

Heavy Duty Protection for Your "Heftiest" Skin Integrity Challenges



Item #	Description
* A8039-29	Geo-Mattress® Atlas, 80″L x 39″W x 7″H
A8042-29	Geo-Mattress® Atlas, 80″L x 42″W x 7″H
A8048-29	Geo-Mattress® Atlas, 80″L x 48″W x 7″H
A8054-29	Geo-Mattress® Atlas, 80″L x 54″W x 7″H
A8060FP-29	Geo-Mattress® Atlas, 80″L x 60″W x 7″H
 F.H. (20//) 750 H. (42// 40// 54// -	1.00//\

Supports 675 lbs. (39"): 750 lbs. (42", 48", 54", and 60").

PressureGuard® Series In-Home Styles Therapy for one. Room for two.



Therapy on Patient Right

Extra width facilitates activities of daily living (ADLs) and care/therapy delivery. Also accommodates couples who want to continue sharing a bed after illness or injury. Fits on our specialty frames or directly on top of standard box springs. When ordering, customer must specify if the dynamic air therapy should be positioned in the **Center** of the mattress (for single user), or on the **Patient's Right Side** (for two users or a single user who cannot/prefers not to move to the center of the mattress). Note: "Patient Right" is with patient in a backlying position.

Note: All models pictured are available in standard hospital width models. Call for model numbers and pricing.



Therapy in Center



Customized lateral rotation 24 hours a day with simple, pendant-controlled made selection and built-in microprocessor-controlled motion system.

PressureGuard Turn Select®

Size/Dimensions	Therapy on Patient Right	Therapy in Center	Weight limit
Twin Long (Dual King*) / 80"L x 38"W		TS8038-29	350 lbs. each user
Full / 75"L x 54"W	TS7554-RT-29	TS7554-CTR-29	350 lbs. each user
Full Long / 80"L x 54"W	TS8054-RT-29	TS8054-CTR-29	350 lbs. each user
Queen / 80"L x 60"W	TS8060-RT-29	TS8060-CTR-29	350 lbs. each user



An easy-access toggle switch for instant changes between alternating pressure and basic lateral rotation.

PressureGuard® APM2 with deluxe control unit

Size/Dimensions	Therapy on Patient Right	Therapy in Center	Weight limit
Twin Long (Dual King*) / 80″L x 38″W		AP28038-29	350 lbs. each user
Full / 75"L x 54"W	AP27554-RT-29	AP27554-CTR-29	350 lbs. each user
Full Long / 80"L x 54"W	AP28054-RT-29	AP28054-CTR-29	350 lbs. each user
Queen / 80"L x 60"W	AP28060-RT-29	AP28060-CTR-29	350 lbs. each user



Powerful results without power. Patented Constant Force Technology uses an air exchange system to provide documented pressure ulcer prevention and treatment without electricity.

PressureGuard CFT®

Size/Dimensions	Therapy across entire surface	Weight limit
Twin Long (Dual King*) / 80"L x 38"W	CF8038-HS-29	500 lbs.
Full / 75"L x 54"W	CF7554-HS-29	750 lbs. total**
Full Long / 80"L x 54"W	CF8054-HS-29	750 lbs. total**
Queen / 80"L x 60"W	CF8060-HS-29	750 lbs. total**
Note: Includes 1 ea. air system renewal hand pump (item #: "SERVICE-00)		

Two Twin Longs side by side would create "Dual King".

Cover cleaning: Lead times:

Non-prorated warranty: 18 months (APM² Deluxe); 2 years (Turn Select), 5 years (CFT). Wipe clean with neutral suds or approved cleaning agents. Avoid bleach. Allow up to 10 working days lead time for all sizes.

All In-Home Style units ship fully assembled in a flat carton via common carrier. **UL Listing:** Easy Air conforms to UL2601-01; APM² Deluxe conforms to UL 544. See individual product brochures for complete technical specifications.

^{** 1} or 2 users.

Here's what's NEW with Transfer Master!

- ✓SUPER Low Bed
- ✓ Single or Dual controls
- ✓Single Massage
- ✓ Floor lighting feature
- ✓220 Volt systems for Export
- ✓Patent Pending split design for easy delivery and assembly



Call 800-475-8122 or visit us at www.transfermaster.com

Customer Testimonials

Here is what a few of our clients have to say about our products and our customer service!



I googled 'full-size electric bed' ...pulled up the Transfer Master site and then everything became very easy. Transfer Master is the epitome of customer service with on-line ease... I dealt with the same rep that was knowledgeable, courteous and extremely responsive to my concerns.. The folks at Transfer Master made what is otherwise usually a very difficult journey ... a pleasure.. *Gail - NH, USA*

.....This bed is truly, outrageously, well made. When I laid my daughter down on it for the first time...... it was the very amused and happy look on her faceShe was and is delighted!.......the fact that this bed raises quite high which saves my back and now I can care for her effortlessly. Not having to stoop even the slightest bit , is a dream.. I have not had a back ache since we got this bed for her. You can bet we will be showing off her new bed and recommending it to everyone we know who could benefit by having one. Thank you so much....! *N. McPherson, PA*